CE xxxx IVD

Precision in Hemostasis

Behnk Elektronik

BE Factor XII Deficient Plasma FXII

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REFERENCIAS

(1) GRIFFIN J.H., COCHRANE C.G.: "Human Factor XII (Hageman factor)" dans « Methods in enzymology », L. Lorand, New York: academic Press, 45, 56-65, 1976.
(2) CAEN J., LARRIEU M-J., SAMMAMA M. : « L'hémostase, méthode d'exploration et diagnostic pratique ». Paris : L'Expansion scientifique, 1975.
(3) IATRIDIS S.G., FERGUSON J.H.: "Effect of physical exercise in blood clotting and fibrinolysis". J. Appl. Physiol., 18, 337-344, 1963.
(4) ANDREW M., PAES B., MILNER R., JOHNSTON M., MITCHELL L., TOLLEFSEN D.M., POWERS P.: "Development of the human coagulation system in the full—term infant". Blood, 70, 1, 165-172, 1987.
(5) REVAK S.D., COCHRANE C.G., GRIFFIN J.H.: "The binding and cleavage characteristics of human Hageman factor during contact activation. A comparison of normal plasma with plasmas deficient in factor XI, prekallikrein, or high molecular weight kininogen". J. Clin.Invest., 59, 1167-1175, 1977.
(6) Young DS. Effects on Clinical Laboratory Tests – Drugs Disease, Herbs & Natural Products, https://clinfx.wiley.com/ aaccweb/aacc/, accessed in September 2021. Published by AACC Press and John Wiley and Sons, Inc.
(7) Secondary Coagulation Standard Lot #5 NIBSC code: SSCLOT5 Instructions for use (Version 3.0, Dated 03/08/2020)
(8) CLSI Document H21-A5: "Collection, transport, and processing of blood specimens for testing plasma-based coagulation assays and molecular hemostasis assays; approved guideline". Fifth edition, 28, 5, 2008.
(9) BACHMANN F.: "Fibrinolysis" in "XIth International congress thrombosis and haemostasis", Thromb. Haemostasis, M. Verstraete, J. Vermylen, H.R. Lijnen, J. Arnout, Louvain: I.S.T.H. and Leuven University Press, 227-265, 1987.
(10) SAMAMA M., CONARD J., HORELLOU M.H., LECOMPTE T.: "Physiologie et exploration de l'hémostase". Paris: Doin, 98-99, 166-173, 1990.
(11) BRANDT J.T., TRIPLETT D.A., ROCK W.A., BOVILL E.G., ARKIR C.F.: "Effect of lupus anticoagulants on the activated partial thromboplastin time". Arch. Pathol. Lab. Med., 115, 109-114, 1991.
(12) SAMPOL J., ARNOUX D., BOUTIERE B.: "Manuel d'hémostase". Paris: Editions scientifiques et médicales Elsevier, 48, 361-362, 1995.
(13) WOODHAMS B., GIRARDOT O., BLANCO M.J., COLESSE G., GOURMELIN Y.: "Stability of coagulation proteins in frozen plasma". Blood Coag. Fibrinolysis, 12, 229-236, 2001.
(14) GIROLAMI A., RUZZON E., LOMBARDI A.M., CABRIO L., RANDI M.L.: "Thrombosis-free surgical procedures in severe (homozygote) factor XII deficiency: report of four additional cases and literature review". Clin. Appl. Thrombosis/Hemostasis, 10, 4, 351-355, 2004.
(15) SCHMAIER A.H., MACCRAE K.R.: "The plasma kallikrein-kinin system: its evolution from contact activation". Journal of Thrombosis and Haemostasis, 5, 2323-2329, 2007.
(16) BLAT Y., SEIFFERT D.: "A renaissance for the contact system in blood coagulation?". Thromb. Haemost., 99, 457-460, 2008

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BIOLABO S.A.S.
Kommanditgesellschaft Behnk Elektronik GmbH & Co.