



The Thrombolyzer XRM is a fully-automated in vitro coagulation analyzer which delivers high quality results and meets the demands of high-volume labs.

## Productivity and Efficiency

Continuous loading and the high on-board capacity of samples, reagents and consumables ensure an uninterrupted workflow and a long walk-away-time. By using two independent samplers for plasma and reagents, combined with an 8-channel measuring system, the analyzer achieves a throughput of 320 tests/h.

## Flexibility and Functionality

Working with the sample rotor ensures maximum achievable flexibility. The sample handling with cap piercing technology for variable primary receptacles provides high user convenience. The positive sample identification, with scanning the barcode directly before sample withdrawal, makes this system absolutely safe.

## Practicality and Reliability

Just a few minutes of daily and weekly maintenance are enough to keep the system up and running. That means lower instrument downtime, less lab personnel expenditure of time and higher processing capacity for patient samples. The XRM, like all Thrombolyzer systems, is well known for reliability with very low operating and maintenance costs.

## Thrombolyzer XRM Specifications:

Sample throughput / with Cap piercing	320 tests/h / 240 tests/h
STAT samples	✓
Time from standby to first PT result	~ 5 min.
Sample identification	Positive, internal barcode reader
Positive reagent identification	X
Clotting, Chromogenic, Immunologic tests	<b>✓</b>
Cap piercing for variable primary receptacles	<b>√</b>
Continuous loading of cuvettes, samples, reagents	<b>✓</b>
Continuous operation / walk-away-time	> 1 hour
Calibration fully automatic	✓
Reaction curve plausibility check	<b>✓</b>
Micro clot detection	✓
Automatic test repetition	✓
Reflex tests	✓
Wavelength	405 nm, 620 nm
Digital light adjustment	✓

Interference detection	$\checkmark$
User-defined tests	✓
Patient data storage	> 100.000
Reaction curve availability	<b>✓</b>
Full reagent and sample result traceability	<b>✓</b>
Reagent positions (cooled) / stirring positions	27(21) / 3
Sample capacity	62
Cuvette capacity	480
Quality control program	✓
Low maintenance effort	Daily 5 min. / weekly 25 min.
Operating system	LINUX
Security system	Three security levels
LIS Interface	Bi-directional
Power consumption	220 VA
Sound level	65 dB (A)
Analyzer dimensions (L x W x H)	106 cm x 60 cm x 78 cm
Weight	58 kg
Mains voltage	100 - 240 V



For further information please feel free to contact us

Phone: +49 (0) 40 - 529861-0 Fax: +49 (0) 40 - 52 98 61-99

info@behnk.de behnk.de Kommanditgesellschaft Behnk Elektronik GmbH & Co. Hans-Böckler-Ring 27 22851 Norderstedt Germany

